


Trade Up Your
UV/Vis



Trade-Up Your UV/Vis Spectrophotometers from Shimadzu or any other vendors for *Better Performance!* Launched in 2019, the UV-1900i UV/Vis Spectrophotometer features several enhancements which will increase productivity and efficiency to your daily routine analyses. **Upgrade NOW and enjoy 30% trade-up discount!**

Email enquiries to sgp@shimadzu.com.sg and quote **2022SGP0601-UV**



UV-1900i UV/Vis Spectrophotometer

High Performance to Meet Diverse Needs:

- Fastest scan speed at 29,000nm/sec, which enables unsurpassed results
- Low stray light of less than 0.02% (220nm) and high photometric repeatability of ± 0.0002 Abs ensures accurate quantitation cross wide concentration range
- LabSolutions UV-Vis file edition is standard. Expansion options with LabSolutions platform for future workstation is possible

Ease of Operation & Network Function through Wireless Connectivity:

- Sleep mode and wakeup function helps reduce warm up time
- Bar-Code Reader and Keyboard ready helps reduce manual entry
- Standalone Ergonomic touch-screen operation or LabSolutions UV-Vis software control

Need higher throughput sample analyses?

HT1500U: Flow-Through Cell Autosampler

Get **20% off** when you purchase with **any** Shimadzu UV/Vis Spectrophotometer

- Flexible configuration to accommodate 4 ml to 40 ml sample tubes
- HT1500U comes with HTA Autosampler Manager, our complementary software for method programming. A simple, easy-to-use keypad allows you to start the analysis by just loading your samples and pushing the START button
- Compatible with Shimadzu UV-1900i, UV-2600i and UV-2700i UV/Vis Spectrophotometers

